




Product Information

Colostrum Preserver

Name:	Colostrum Preserver	
Pack Sizes/Codes:	5L / DBCP5 20L /.DBC20 200L / DBCP200	
Main Points:	<ul style="list-style-type: none"> Formulated to inhibit bacterial growth and prolong the useable life of all milk based products used for calf rearing Will preserve colostrum and calf milk for 3 -4 weeks 	
Status:	OTC	
Ingredients:	<ul style="list-style-type: none"> Methanoic Acid 10-15% (feed grade) 	
Indications:	<ul style="list-style-type: none"> Colostrum Preserver will preserve colostrum and calf milk for three to four weeks by slowing its normal degradation, preventing foul odours and loss of nutritional value. The product can be effectively used for the prolonged preservation of milk replacers, colostrum, whole milk, skim milk, whey and other milk products fed to animals. By using AHD Colostrum Preserver, the fat, protein and lactose quality of the milk products are preserved and their nutritional value is maintained. 	
Action:	Inhibits bacterial growth and prolongs the useable life of milk based products	
Administration:		
Dosage:	<ul style="list-style-type: none"> Mix AHD Colostrum Preserver at the rate of 200ml to 20 litres fresh colostrum milk. For best results: Add Colostrum Preserver on a regular basis as fresh milk product is added to the mix. Separation will occur over the holding period, this is a natural process, regular agitation will assist with the preservation process. NOTE: AHD Colostrum Preserver treated milk products should NOT be heated prior to feeding to animals. 	
Special Precautions:	<p>AHD Colostrum milk should NOT be heated prior to feeding to animals.</p> <p>This product causes burns to eyes, skin and mucous membranes. May be corrosive to metals.</p> <p>Take precautions to avoid all contact; wear protective goggles, impervious gloves, apron and boots.</p> <p>Work in well-ventilated location</p>	
Storage:	Keep out of reach of children.	
Registration:	Exempt	
Category:	Dietary supplement	



Animal Health Direct Limited

PO Box 8015 | Havelock North | 4157 | t: 06 873 3611 | f: 06 873 3616
 e: sales@ahdltd.co.nz | ahdltd.co.nz

