F-Grip Plastic Injector Range

- ¬ Lightweight, accurate and easy to use.
- ¬ Superior comfort grip and easy-off needle nut.
- \neg Supplied with 2 needles.
- ¬ Sterilisable by boiling.
- ¬ No spare parts available.




MODEL	DOSE RANGE	INCREMENTS	CODE
1ml	0.1ml – 1ml	0.1ml	PAS817
3ml	0.25ml – 3ml	0.25ml	PAS818
5ml	0.25ml – 5ml	0.25ml	PAS819
10ml	0.5ml – 10ml	0.5ml	PAS820

SPECIAL NOTE

In North America these applicators are supplied without needles, but with a feed tube and a universal needle nut to suit aluminium, plastic and brass hub needles.

Codes:	PAS817F – 1ml
	PAS819F – 5ml
	PAS818F – 3ml
	PAS820F - 10ml



Customise this instrument with your brand colour and logo. Contact us for details.

Njphillips & com Our expertise is in your hands.

LOCKED BAG 8, GOSFORD NSW, AUSTRALIA 2250 P +61 2 4340 2044 F +61 2 4340 1991 | E njp1@njphillips.com.au

The information contained in this document is legally privileged and confidential intended only for the use of the individual or entity. If the receiver of this document is not the intended recipient, the receiver is hereby notified that any dissemination, distribution or copy of this information is strictly prohibited. If this information is received in error, please notify us immediately. As our policy is one of continuous improvement, the manufacturer reserves the right to alter these specifications at any time. All products produced by NJ Philips Pty Limited are identified by a unique batch number which is affixed to the product to allow traceability by the manufacturer and must not be removed if product integrity is to be maintained NoI all part numbers may be available.